

GM1361 Temperature & Humidity Meter



[Click here to view product](#)

Product Code	XE-4012-001
Model	GM1361 Temperature & Humidity Meter
General Description	Supplied in hard carry case
Type	K-type thermocouple temperature measuring
Accuracy of temperature	±1°C or ±1.8°C
Accuracy of humidity	± 3%(30-95%); ± 5%(10-30%)
Battery	9V Battery
Dimensions	160*70*30mm
Function #1	Large LCD temperature & humidity dual display
Function #2	Data hold

Function #3	Low battery indication
Function #4	Auto /Manual power shut off
Function #5	Max/ Min function
Input Range - Type K	-10°C ~ 50°C (-30°C ~ 1000 °C K-type)
Humidity range	5%RH ~ 98%RH
Unit selectable	°C / °F selection
Sampling rate	2.5times/sec