

GM1361 Temperature & Humidity Meter



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Currently awaiting new stock - pre-order only.

You will be notified when your order is ready to ship.

Specifications

Product Code	XE-4012-001
Model	GM1361 Temperature & Humidity Meter
General Description	Supplied in hard carry case
Type	K-type thermocouple temperature measuring
Accuracy of temperature	$\pm 1^{\circ}\text{C}$ or $\pm 1.8^{\circ}\text{C}$
Accuracy of humidity	$\pm 3\%$ (30-95%); $\pm 5\%$ (10-30%)
Battery	9V Battery
Dimensions	160*70*30mm
Function #1	Large LCD temperature & humidity dual display
Function #2	Data hold
Function #3	Low battery indication
Function #4	Auto /Manual power shut off
Function #5	Max/ Min function
Input Range - Type K	$-10^{\circ}\text{C} \sim 50^{\circ}\text{C}$ ($-30^{\circ}\text{C} \sim 1000^{\circ}\text{C}$ K-type)
Humidity range	5%RH ~ 98%RH
Unit selectable	$^{\circ}\text{C}$ / $^{\circ}\text{F}$ selection
Sampling rate	2.5times/sec