

Fixed Racing Tyre Temperature Probe



[Click here to view product](#)

Product Code	XE-3891-001
General Description	Can be used with any digital thermometer with a type k mini input socket
Sensor Type	Fast response fixed racing tyre probe with a max depth of six millimetres
Thermocouple Type	K
Calibration	IEC
Probe Length	max depth of 6mm
Cable Length	1m
Cable Termination	Miniature Plug
Max. Temperature	+250°C
Min. Temperature	-50°C

Features #1	Junction Type: Grounded
-------------	-------------------------

