



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

PFA Insulated Flat Pair IEC Exposed Junction Thermocouple with Miniature Plug - Type K - 1/0.3mm Strand Size - 1M Length



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Good mechanical strength and flexibility, resistant to oils, acids and other adverse fluids.

Need Your Temperature Sensor calibrated?

Traceable to ISO17025 (UKAS) calibration standards, Labfacility offer a choice of a 3 or 5 point traceable calibration using our in-house calibration facility, you can select any temperature points between -10°C & +200°C.

Please be aware that if you add a calibration option to your order you will be contacted regarding delivery information

Specifications

Specifications

Product Code	XE-3484-001
Sensor Type	Fast response, welded exposed junction
Thermocouple Type	K IEC
Cable Connection	Terminated with a miniature flat pin plug
Cable Length	1M
Core / Strands	1/0.3mm
Termination Type	Miniature Plug
Standards Met	(cores & jacket) to IEC-584-3
Insulation Material	Insulation rating: -75°C to 260°C
Series	Exposed Junctions