

PFA Insulated Flat Pair IEC Exposed Junction Thermocouple with Miniature Plug - Type K - 1/0.3mm Strand Size - 1M Length



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.



Good mechanical strength and flexibility, resistant to oils, acids and other adverse fluids.

Specifications

Specifications

Fast response, welded exposed junction
. det respected, traided apparatus
K IEC
Terminated with a miniature flat pin plug
1M
1/0.3mm
Miniature Plug
(cores & jacket) to IEC-584-3
Insulation rating: -75°C to 260°C
Exposed Junctions