



# **LABFACILITY**

TEMPERATURE & PROCESS TECHNOLOGY

## **PFA Insulated Flat Pair IEC Exposed Junction Thermocouple with Miniature Plug - Type K - 1/0.3mm Strand Size - 3M Length**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Good mechanical strength and flexibility, resistant to oils, acids and other adverse fluids.

### **Need Your Temperature Sensor calibrated?**

Traceable to ISO17025 (UKAS) calibration standards, Labfacility offer a choice of a 3 or 5 point traceable calibration using our in-house calibration facility, you can select any temperature points between -10°C & +200°C.

**Please be aware that if you add a calibration option to your order you will be contacted regarding delivery information**

### **Specifications**

## **Specifications**

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | XE-3486-001                               |
| <b>Sensor Type</b>         | Fast response, welded exposed junction    |
| <b>Thermocouple Type</b>   | K IEC                                     |
| <b>Cable Connection</b>    | Terminated with a miniature flat pin plug |
| <b>Cable Length</b>        | 3M  |
| <b>Core / Strands</b>      | 1/0.3mm                                   |
| <b>Termination Type</b>    | Miniature Plug                            |
| <b>Standards Met</b>       | (cores & jacket) to IEC-584-3             |
| <b>Insulation Material</b> | Insulation rating: -75°C to 260°C         |
| <b>Series</b>              | Exposed Junctions                         |