


## **PFA Insulated Flat Pair IEC Exposed Junction Thermocouple with Miniature Plug - Type K - 1/0.5mm Strand Size - 1M Length**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Good mechanical strength and flexibility, resistant to oils, acids and other adverse fluids.

### Specifications

## Specifications

<b>Product Code</b>	XE-3488-001
<b>Sensor Type</b>	Fast response, welded exposed junction
<b>Thermocouple Type</b>	K IEC
<b>Cable Connection</b>	Terminated with a miniature flat pin plug
<b>Cable Length</b>	1M
<b>Core / Strands</b>	1/0.3mm
<b>Termination Type</b>	Miniature Plug
<b>Standards Met</b>	(cores & jacket) to IEC-584-3
<b>Insulation Material</b>	Insulation rating: -75°C to 260°C