


PFA Insulated Flat Pair IEC Exposed Junction Thermocouple with Miniature Plug - Type K - 1/0.5mm Strand Size - 5M Length



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Good mechanical strength and flexibility, resistant to oils, acids and other adverse fluids.

Specifications

Specifications

Product Code	XE-3491-001
Sensor Type	Fast response, welded exposed junction
Thermocouple Type	K IEC
Cable Connection	Terminated with a miniature flat pin plug
Cable Length	5M
Core / Strands	1/0.3mm
Termination Type	Miniature Plug
Standards Met	(cores & jacket) to IEC-584-3
Insulation Material	Insulation rating: -75°C to 260°C
Series	Exposed Junctions