

Terminal Blocks - Type J IEC



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Terminal blocks come in many different types and sizes to suit your terminal head.

Check carefully the dimensions of the fixing holes on your terminal head before choosing the correct block.

Way

Order Code

Colour

Diameter

Fixings

Fits Head

2

XE-3827-001

Black

42mm

33mm

KNE, KNP, KPP, B, SCH4

4

XE-3828-001

Black

42mm

33mm

KNE, KNP, KPP, B, SCH4

Specifications

Material	Plastic
Fixing (Centres)	33mm
Colour	Black
Dimensions	42mm Diameter
To Fit Head	KNE, KNP, KPP, B, SCH4
Sizes	2 or 4 Way
Thermocouple Type	J