

KNP Hinged Lid Terminal Head



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Many alternative types of terminal heads are available to meet the requirements of various applications.

Variations exist in size, material, accommodation, resistance to media, resistance to fire or even explosion and in other parameters.

The KNP Hinged Lid Terminal Head is a silver-painted aluminium enclosure designed to protect electrical connections, particularly for temperature sensors and transmitters in industrial environments. Its hinged cap design provides easy access for maintenance and installation.

Specifications

Specifications

Head Type	KNP Hinged Lid Terminal Head
Material	Aluminium Silver Painted
Cable Entry	M20 x 1.5mm
Fixing (Centres)	33mm & 40mm
Max. Temperature	170°C
IP Rating	IP67 rated
Process Entry	1/2" BSPP and 3/4" BSPP
Cap Retaining Type	Hinged Cap