

## KF Wall Mounted Terminal Head



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Many alternative types of terminal heads are available to meet the requirements of various applications.

Variations exist in size, material, accommodation, resistance to media, resistance to fire or even explosion and in other parameters.

The KF Wall Mounted Terminal Head is an industrial terminal head made from die-cast aluminium with a black or silver painted finish, designed to provide robust protection for sensing elements in temperature measurement applications.

### Specifications

## Specifications

Head Type	KF Wall Mounted Terminal Head
Material	Aluminium Silver Painted
Cable Entry	M20 x 1.5mm
Fixing (Centres)	33mm & 40mm
Max. Temperature	170°C
IP Rating	IP67 rated
Process Entry	1/2" BSPP
Cap Retaining Type	Screw Cap and retaining chain