

Lascar EL-MOTE-TP - Wireless Cloud-Connected Temperature Data Logger with Thermistor Probe



[Click here to view product](#)

Product Code	XE-4201-001
General Description	Perfect for both compact systems with just a few sensor requirements to corporate solutions with thousands of devices around the globe
Type	Wireless Cloud-Connected Temperature Data Logger with Thermistor Probe
Communications Protocol	Wireless
Accuracy	$\pm 0.8^{\circ}\text{C}$ $\pm 1.5^{\circ}\text{F}$
Dimensions	93 x 93 x 32mm 3.7 x 3.7 x 1.3"
Length	Probe Cable Length 3m (118.1")
Input / Ranges	-40 to +125°C -40 to +257°F
Power Consumption	Up to 2 years 4 x 1.5V AA cells

Accessories	Mains adaptor for permanent power EL-MOTE-PSU
-------------	---