

## MA Terminal Head



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Many alternative types of terminal heads are available to meet the requirements of various applications.

Variations exist in size, material, accommodation, resistance to media, resistance to fire or even explosion and in other parameters.

The MA Terminal Head is a compact, robust, micro-style connection head primarily used with industrial temperature sensors in applications where space is limited, aluminium silver painted.

### Specifications

## Specifications

Product Code	XE-4129-001
Head Type	MA Terminal Head
Material	Aluminium
Cable Entry	M16 x 1.5mm
Fixing (Centres)	20mm
Max. Temperature	170°C
IP Rating	IP66 rated
Process Entry	1/8" BSPP
Cap Retaining Type	Cap and retaining screw