

## N Terminal Head



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Many alternative types of terminal heads are available to meet the requirements of various applications.

Variations exist in size, material, accommodation, resistance to media, resistance to fire or even explosion and in other parameters.

## Specifications

<b>Product Code</b>	XE-4130-001
<b>Head Type</b>	N Terminal Head
<b>Material</b>	Nylon Black Head
<b>Cable Entry</b>	M20
<b>Fixing (Centres)</b>	33mm
<b>Max. Temperature</b>	90°C
<b>IP Rating</b>	IP54 rated
<b>Process Entry</b>	1/8" BSPP
<b>Cap Retaining Type</b>	Screw Cap Lid with Rubber Gasket