

Pt1000 Class B Sensors with Teflon insulated lead - Type PRT - 4.0mm Diameter - 75mm Length



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-3683-001
General Description	Insulated element in a 4.0mm diameter 316 stainless steel probe
Sensor Type	Pt1000 (1000 Ohms @ 0°C) to IEC 751, Class B
Element	Pt1000
Probe Diameter	4.0mm
Probe Length	75mm
Sheath Material	316 stainless steel probe
Cable Length	1M
Cable Type	Teflon® insulated twisted leads
Cable Termination	Bare tails
Max. Temperature	200°C
Min. Temperature	-50°C
Number of Wires	2
Series	Pt1000 Sensors