



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## **Pt1000 Class B Sensors with Teflon insulated lead - Type PRT - 4.0mm Diameter - 125mm Length**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Specifications

**Specifications**

<b>Product Code</b>	XE-3684-001
<b>General Description</b>	Insulated element in a 4.0mm diameter 316 stainless steel probe
<b>Sensor Type</b>	Pt1000 (1000 Ohms @ 0°C) to IEC 751, Class B
<b>Element</b>	Pt1000
<b>Probe Diameter</b>	4.0mm
<b>Probe Length</b>	125mm
<b>Sheath Material</b>	316 stainless steel probe
<b>Cable Length</b>	1M
<b>Cable Type</b>	Teflon® insulated twisted leads
<b>Cable Termination</b>	Bare tails
<b>Max. Temperature</b>	200°C
<b>Min. Temperature</b>	-50°C
<b>Number of Wires</b>	2
<b>Series</b>	Pt1000 Sensors