

## Legionella Thermometer with Integral Water Temp Timer



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

MM2008 Single input Thermometer with an integral 1 minute & 2 minute timer.

This Thermometer is primarily developed for use in Legionella risk management and offers reassurance that the correct temperature is taken when monitoring hot and cold water temperatures.

### Features

Single Input Thermocouple

Integral Timer for 1 and 2 minute countdown

Measurement Range: -200 to +1,372°C (probe type dependant)

Accuracy: -  $\pm 0.1\%$  of Reading  $\pm 0.2^\circ\text{C}$

C/F Selectable

IP67 Waterproof Casing

Low Battery Indication

Selectable Thermocouple Types: K & T

Self-Calibrating

Rubber Boot included for extra protection against damage

Includes 1 x PP3 Battery

Requires: A Thermocouple Probe/Sensor

### Specifications

## Specifications

|                     |   |
|---------------------|---|
| Product Code        | XE-4091-001   |
| General Description | Measurement Range: -200 to +1,372°C (probe type dependant)          |
| Type                | Single input Thermometer with an integral 1 minute & 2 minute timer |
| Accuracy            | $\pm 0.1\%$ of Reading $\pm 0.2^\circ\text{C}$                      |
| Readings            | C/F Selectable  |
| Power Consumption   | Low Battery Indication  |
| Power Supply        | Includes 1 x PP3 Battery  |



|                          |  |
|--------------------------|--|
| <b>Standards Met</b>     | IP67 Waterproof Casing                                   |
| <b>Thermocouple Type</b> | K & T  |
| <b>Thermometer</b>       | Self-Calibrating   |
| <b>Termination Type</b>  | Requires a Thermocouple Probe/Sensor                     |
| <b>Accessories</b>       | Rubber Boot included for extra protection against damage |