

Legionella Thermometer with Integral Water Temp Timer



[Click here to view product](#)

| | |
|----------------------------|---|
| Product Code | XE-4091-001 |
| General Description | Measurement Range: -200 to +1,372°C (probe type dependant) |
| Type | Single input Thermometer with an integral 1 minute & 2 minute timer |
| Accuracy | ± 0.1% of Reading ± 0.2°C |
| Readings | C/F Selectable |
| Power Consumption | Low Battery Indication |
| Power Supply | Includes 1 x PP3 Battery |
| Standards Met | IP67 Waterproof Casing |
| Thermocouple Type | K & T |
| Thermometer | Self-Calibrating |

| | |
|-------------------------|--|
| Termination Type | Requires a Thermocouple Probe/Sensor |
| Accessories | Rubber Boot included for extra protection against damage |

