

Legionella Thermometer with Integral Water Temp Timer



[Click here to view product](#)

Product Code	XE-4091-001
General Description	Measurement Range: -200 to +1,372°C (probe type dependant)
Type	Single input Thermometer with an integral 1 minute & 2 minute timer
Accuracy	± 0.1% of Reading ± 0.2°C
Readings	C/F Selectable
Power Consumption	Low Battery Indication
Power Supply	Includes 1 x PP3 Battery
Standards Met	IP67 Waterproof Casing
Thermocouple Type	K & T
Thermometer	Self-Calibrating

Termination Type	Requires a Thermocouple Probe/Sensor
Accessories	Rubber Boot included for extra protection against damage