

**Brass Flange - 50mm Diameter - 1/2" BSPT
Thread - 40mm PCD**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Brass Flanges are crucial components in piping systems, offering robust and reliable connections for various industrial applications.

Their corrosion resistance, durability, and strength make them suitable for demanding environments in industries such as oil and gas, chemical processing, power generation, food and beverage, and pharmaceuticals.

Proper selection, installation, and maintenance of stainless steel flanges ensure efficient and long-lasting performance in any piping system.

Specifications

Specifications

Product Code	XE-5452-001
General Description	Flanges can be used to secure the sensor assembly into the process.
Material	Brass
Dimensions	50mm dia, 40mm PCD
Sizes	1/2" BSPT