

RTD Thermometer (dual input) Model CHY 804U



[Click here to view product](#)

Product Code	XE-3932-BUNDLE
Type	RTD Thermometer (dual input)
Accuracy	$\pm(0.05\%rdg + 0.2^{\circ}C)$ on $^{\circ}C$ scale $\pm(0.05\%rdg + 0.4^{\circ}F)$ on $^{\circ}F$ scale
Resolution	0.1 $^{\circ}C$ or 0.2 $^{\circ}F$
Logging Interval	1 times/second
Input / Ranges	Pt385(100ohms): -200 $^{\circ}C$ to 800 $^{\circ}C$, (-328 $^{\circ}F$ to 1472 $^{\circ}F$) Pt3916/ Pt3926(100ohms): -200 $^{\circ}C$ to 630 $^{\circ}C$ (-328 $^{\circ}F$ to 1166 $^{\circ}F$)
Operating Temperature	0 $^{\circ}C$ to 50 $^{\circ}C$ at<70% R.H
Weight	approx. 255g including batteries.
Power Supply	1.5V x 4pcs AAA size.