


## RTD Temperature Sensor with Hirschmann Connector - 100mm Probe Length



**'HIRSCHMANN PT100 PROBE ASSEMBLY'**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

A cost effective 6mm diameter Pt100 sensor assembly housed in a 316 Stainless Steel body with 1/2" BSP process connection and fitted din connector with connection lead.

### Specifications

## Specifications

Product Code	XE-5683-001
General Description	RTD Temperature Sensor & Hirschmann Connector
Sensor Type	PT100 Class A
Process Connection	1/2" BSP
Probe Diameter	6.0mm
Probe Length	100mm
Cable Length	1m
Cable Type	PVC Insulated
Max. Temperature	250°C
Min. Temperature	-50°C
IP Rating	IP65
Connection Type	DIN connector Screw terminal
Number of Wires	4 wire connection
Series	Hirschmann Temperature Sensor