


## **IEC Type K 0.7kg Pull Button Magnet Thermocouple, PTFE Insulated with Type K Miniature Plug - 1M Length**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Specifications****Specifications**

<b>Product Code</b>	XE-3460-001
<b>General Description</b>	Ideal for ferrous metal surface measurement
<b>Thermocouple Type</b>	K
<b>Thermocouple Junction</b>	Grounded at tip
<b>Magnet</b>	1/2" diameter button type
<b>Cable Length</b>	1 or 2 metres
<b>Cable Type</b>	PTFE insulated twin-twist wires
<b>Cable Termination</b>	Miniature or Standard Fitted Plug
<b>Max. Temperature</b>	+250°C (plug 200°C)
<b>Series</b>	Magnet Thermocouple