

IEC Type K 0.7kg Pull Button Magnet Thermocouple, PTFE Insulated with Type K Miniature Plug - 2M Length



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-5617-001
General Description	Ideal for ferrous metal surface measurement
Thermocouple Type	K
Thermocouple Junction	Grounded at tip
Magnet	1/2" diameter button type
Cable Length	1 or 2 metres
Cable Type	PTFE insulated twin-twist wires
Cable Termination	Miniature or Standard Fitted Plug
Max. Temperature	+250°C (plug 200°C)
Series	Magnet Thermocouple