

IEC Type K 3kg Pull Magnet Thermocouple, PFA Insulated Cable with Stainless Steel Over-Braid Terminating in Miniature Plug - 2M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-3470-001
General Description	Rapid, reliable temperature sensing on ferrous surfaces
Thermocouple Type	K
Thermocouple Junction	Grounded measuring junction
Magnet	1" diameter (3kg pull)
Cable Length	2M
Cable Type	stainless steel over-braided PFA Teflon®
Cable Termination	Miniature Plug
Max. Temperature	+250°C (200°C Plug)
Standards Met	IEC
Series	Magnet Thermocouple