


## ANSI 0.25mm Diameter Fast Response Mineral Insulated Thermocouple with Threaded Pot & PFA ANSI Tails - Type K Grounded Junction - 150mm Length



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.



custom sensor form

Please use our s to submit your requirements for various lengths of MI and Cable lengths.

## **Need Your Temperature Sensor calibrated?**

Traceable to ISO17025 (UKAS) calibration standards, Labfacility offer a choice of a 3 or 5 point traceable calibration using our in-house calibration facility.

Please be aware that if you add a calibration option to your order you will be contacted regarding delivery information.

## **Specifications**

## **Specifications**

XE-0816-001
Highly flexible, sheath can be bent/formed to suit many applications and
processes Type K ANSI MD STYLE probes 0.25mm
Single element, grounded junction
150mm
Stainless steel 316 sheath
M8 x 1.0mm fine pitch threaded pot seal (200°C)
100mm
PFA insulated twin twisted ANSI
7/0.2mm
260°C
-40°C
A 1%+ Length Tolerance applies
ANSI MC96.1