



## PCB Mounting Thermocouple Connector Socket FMTC-CU-PCB Type Cu



**Miniature PCB Socket**  
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**PCB Mounting Thermocouple Connector Socket FMTC-CU-PCB Type CU**

Specially designed for PCBs, solid pins for direct board mounting & provision for CJC sensor

**Specifications****Specifications**

<b>Product Code</b>	XE-1379-001
<b>General Description</b>	Versatile, high quality circuit board mounting miniature socket
<b>Thermocouple Type</b>	IEC, ANSI, JIS Type CU
<b>Connector Type</b>	PCB socket
<b>Connector Size</b>	Miniature
<b>Colour</b>	White
<b>Max. Temperature</b>	220°C
<b>Pins</b>	Brass
<b>+Positive Contact</b>	Copper
<b>-Negative Contact</b>	Copper
<b>Standards Met</b>	BS EN 50212:1997. RoHS
<b>Series</b>	Miniature Connector