

0.75mm Diameter with Bare Tails ? Insulated or Grounded



[Click here to view product](#)

Product Code	XE-0191
Sensor Type	Mineral Insulated
Thermocouple Type	K
Element/hot junction	Single element, junction insulated or grounded from sheath
Probe Diameter	0.75mm
Probe Length	100,150, 250mm
Sheath Material	K - Stainless steel 316
Termination Type	Fortfix sealed with = +VE leg 15mm long and -VE leg 10mm long
Max. Temperature	+650°C
Min. Temperature	-40°C

Tolerance	Class 1
Standards Met	In accordance with IEC 584