

Digital Cooking Thermometer Timer With Alarm & rear magnet



[Click here to view product](#)

Product Code	XE-4350-001
Model	LDT-100
General Description	Digital Cooking Thermometer Timer With Alarm & rear magnet
Type	Timer-Thermometer 2 in 1
Accuracy	+/-1°C (0-10°C) /-1.5°C(101-25°C)
Readings	Timer range: 99 Min 59 Sec
Resolution	0.1°C(0.2°F)
High Alarm Level	High temperature alarm function
Battery	1 x 1.5V(AAA)
Cable Length	Wire: 120cm length

Length	10m
Colour	Metal finish
Dimensions	6.5 x 7 x 1.7cm
Functions	°C / °F switch
Temperature Range	0-250°C / 32-482°F
Probe Material	304 stainless steel
Weight	100g
Construction	Eco-friendly ABS plastic + Stainless steel front shell
Important Information	Hanging & flip stand