




# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## **IEC Type K 1.9kg Pull Button Magnet Thermocouple, PFA Insulated Cable with Stainless Steel Over-Braid Terminating in Miniature Plug - 1M**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Specifications****Specifications**

<b>Product Code</b>	XE-0079-001
<b>General Description</b>	1.9kg Pull Button Magnet Thermocouple
<b>Sensor Type</b>	Magnet Thermocouple
<b>Thermocouple Type</b>	K
<b>Magnet</b>	19.1mm Diameter / 12.7mm height. Max Pull: 1.9kg
<b>Cable Length</b>	1M
<b>Cable Type</b>	PFA Flat Pair with stainless steel overbraid
<b>Core / Strands</b>	7/0.2mm
<b>Cable Termination</b>	Miniature Plug
<b>Max. Temperature</b>	250°C
<b>Min. Temperature</b>	-50°C
<b>Standards Met</b>	IEC