




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

IEC Type K 1.9kg Pull Button Magnet Thermocouple, PFA Insulated Cable with Stainless Steel Over-Braid Terminating in Miniature Plug - 2M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-0080-001
General Description	1.9kg Pull Button Magnet Thermocouple
Sensor Type	Magnet Thermocouple
Thermocouple Type	K
Magnet	19.1mm Diameter / 12.7mm height. Max Pull: 1.9kg
Cable Length	2M
Cable Type	PFA Flat Pair with stainless steel overbraid
Core / Strands	7/0.2mm
Termination Type	Miniature Plug
Max. Temperature	250°C
Min. Temperature	-50°C
Standards Met	IEC