

## **ANSI Silicone Rubber Patch Thermocouple, PFA Twin Twisted - Type K With Miniature Plug Termination - 1M**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Specifications****Specifications**

|                            |                                                               |
|----------------------------|---------------------------------------------------------------|
| <b>Product Code</b>        | XE-0573-001                                                   |
| <b>General Description</b> | Resistant to many oils & chemicals                            |
| <b>Sensor Type</b>         | Thermocouple junction encapsulated in a silicone rubber patch |
| <b>Thermocouple Type</b>   | K                                                             |
| <b>Cable Length</b>        | 1M                                                            |
| <b>Cable Type</b>          | PFA insulated twin-twist wire                                 |
| <b>Core / Strands</b>      | 1/0.2mm                                                       |
| <b>Termination Type</b>    | Miniature Plug                                                |
| <b>Max. Temperature</b>    | 150°C                                                         |
| <b>Min. Temperature</b>    | -150°C                                                        |
| <b>Standards Met</b>       | ANSI MC96.1                                                   |
| <b>Colour</b>              | Positive: Red / Negative: Yellow                              |
| <b>Dimensions</b>          | Length 40mm x Width 13mm x Height 5mm Rubber patch            |
| <b>Series</b>              | Rubber Patch Thermocouple                                     |