

## **ANSI Type K 0.7kg Pull Button Magnet Thermocouple, PFA Insulated with Miniature Plug - 0.5M**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Specifications****Specifications**

|                              |                                             |
|------------------------------|---------------------------------------------|
| <b>Product Code</b>          | XE-0581-001                                 |
| <b>General Description</b>   | Ideal for ferrous metal surface measurement |
| <b>Thermocouple Type</b>     | K                                           |
| <b>Thermocouple Junction</b> | Grounded at tip                             |
| <b>Magnet</b>                | 1/2" diameter button type, 0.7kg Pull       |
| <b>Cable Length</b>          | 0.5M                                        |
| <b>Cable Type</b>            | PFA insulated twin-twist wires              |
| <b>Cable Termination</b>     | Miniature Fitted Plug                       |
| <b>Max. Temperature</b>      | +250°C (plug 200°C)                         |
| <b>Standards Met</b>         | ANSI MC96.1                                 |
| <b>Colour</b>                | + Yellow / - Red                            |
| <b>Series</b>                | Magnet Thermocouple                         |