



# **LABFACILITY**

TEMPERATURE & PROCESS TECHNOLOGY

## **ANSI Bolt Thermocouple, Glassfibre Stainless Steel Over Braided Cable - Type K - M10 Thread - 13mm Bolt**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Specifications****Specifications**

<b>Product Code</b>	XE-0586-001
<b>Thermocouple Type</b>	K
<b>Thermocouple Junction</b>	Bolts have grounded junctions
<b>Cable Length</b>	2M
<b>Cable Type</b>	Glassfibre stainless steel over braided cable
<b>Cable Termination</b>	Bare tails
<b>Max. Temperature</b>	250°C
<b>Standards Met</b>	ANSI MC96.1
<b>Bolt thread/pitch</b>	M10 x 1mm
<b>Colour</b>	Positive: Red / Negative: Yellow
<b>Length</b>	13mm Bolt
<b>Series</b>	Bolt Thermocouple