



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

JIS Type K 1.9kg Pull Button Magnet Thermocouple, PFA Insulated Cable with Stainless Steel Over-Braid Terminating in Miniature Plug - 2m



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-0625-001
General Description	1.9kg Pull Button Magnet Thermocouple
Thermocouple Type	K JIS
Magnet	19.1mm Diameter / 12.7mm height - Max pull: 1.9 KG
Cable Length	2m
Cable Type	PFA Insulated Cable with stainless steel overbraid 7/0.2mm
Termination Type	Miniature Plug
Max. Temperature	250°C
Min. Temperature	-50°C
Standards Met	JIS (Japanese Industrial Standards)