



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

ANSI Fabricated Grounded Thermocouple in Stainless Steel Tube, Glassfibre Stainless Steel Overbraided Cable - Type K - 3.18mm Diameter - 25mm Length



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-0633-001
General Description	316 stainless steel tube in various sizes
Thermocouple Type	K
Thermocouple Junction	Grounded at tip for fast response
Cable Length	2M
Cable Type	Glassfibre stainless steel overbraided cable
Cable Termination	Bare tails
Max. Temperature	350°C
Min. Temperature	-60°C
Standards Met	ANSI MC96.1
Series	Fabricated Tube Thermocouple