

**ANSI Type K 4.5kg Pull Power (Horseshoe)
Magnet Thermocouple, PFA Insulated Cable with
Stainless Steel Over-Braid Terminating in
Standard Plug - 1M**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-0668-001
General Description	4.5kg Pull Power (Horseshoe) Magnet Thermocouple
Sensor Type	Magnet Thermocouple
Thermocouple Type	K
Magnet	30 x 20 x 20mm (L x H x W). Max Pull: 4.5kg
Cable Length	1M
Cable Type	PFA Teflon® insulated with stainless steel over-braid
Core / Strands	7/0.2mm
Termination Type	Standard Plug
Max. Temperature	250°C
Min. Temperature	-50°C
Standards Met	ANSI MC96.1
Series	Magnet Thermocouple