

# JIS Type K 4.5kg Pull Power (Horseshoe) Magnet Thermocouple, PFA Insulated Cable with Stainless Steel Over-Braid Terminating in Miniature Plug - 2m



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Specifications****Specifications**

<b>Product Code</b>	XE-0675-001
<b>General Description</b>	4.5kg Pull Power (Horseshoe) Magnet Thermocouple
<b>Sensor Type</b>	K JIS
<b>Magnet</b>	30 x 20 x 20mm (L x H x W). Max Pull: 4.5kg
<b>Cable Length</b>	2m
<b>Cable Type</b>	Glassfibre insulated lead with stainless steel overbraid
<b>Core / Strands</b>	7/0.2mm
<b>Termination Type</b>	Miniature Plug
<b>Max. Temperature</b>	+250°C
<b>Min. Temperature</b>	-50°C
<b>Standards Met</b>	JIS (Japanese Industrial Standards)