


# JIS Type K General Purpose Thermocouple Probe, Glassfibre insulated Cable with Stainless Steel Over-Braid Terminating in Bare tails, Miniature or Standard Plug



[Click here to view product](#)

Product Code	XE-0768
General Description	Rugged thermocouple with stainless steel sheath
Thermocouple Type	K
Thermocouple Junction	Grounded hot junction for fast response
Probe Diameter	4.5mm
Probe Length	150mm
Cable Length	1 or 2m
Cable Type	Glassfibre insulated lead with stainless steel overbraid
Core / Strands	7/0.2mm
Termination Type	Bare tails, Miniature or Standard Plug (with cable clamps)

<b>Max. Temperature</b>	+350°C
<b>Min. Temperature</b>	-60°C
<b>Standards Met</b>	JIS