

ANSI Nozzle Thermocouple, Glassfibre insulated Cable with Stainless Steel Over-Braid - Type K Terminating in Bare Tails - 1M Length



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.



Specifications

Specifications

XE-0869-001
Thermocouple housed in a stainless steel threaded bolt
K
Grounded hot junction for fast response, M6x1mm thread
1M
Glassfibre insulated lead with stainless steel overbraid
7/0.2mm
Bare tails
+350°C
-60°C
ANSI MC96.1
Nozzle Thermocouple