

**Washer Thermocouple, Glassfibre insulated
Cable with Stainless Steel Over-Braid - Type K JIS
with Bare tails Termination and 6.0mm Washer
Hole Size**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-0976-001
General Description	Thermocouple housed in a washer eyelet
Thermocouple Type	K JIS
Thermocouple Junction	Grounded tip for quick response
Cable Length	2m
Cable Type	Glassfibre insulated, stainless steel overbraided
Core / Strands	7/02.mm
Termination Type	Bare Tails
Max. Temperature	+350°C
Standards Met	JIS (Japanese Industrial Standards)
Washer	6.0mm clearance hole