

IEC Fine Wire Disc Thermocouple, PTFE with Miniature Plug - Type J - 5M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

Specifications

Product Code	XE-3760-001
General Description	Exposed junction thermocouple with Copper Disc brazed at tip
Sensor Type	Copper disc
Thermocouple Type	J
Thermocouple Junction	Grounded measuring junction
Cable Length	5M
Cable Type	1/0.2mm PTFE insulated twisted wire
Termination Type	Miniature flat pin plug, colour coded in accordance IEC 584
Standards Met	IEC