

Flat Tip Probe ? Type RTD / PRT



[Click here to view product](#)

Product Code	XE-3697
General Description	316 stainless steel probe, 4.0mm diameter (6mm max. across flat area of tip) x 150mm long overall
Sensor Type	Pt100 (100 Ohms @ 0°C) to IEC 751, Class B
Probe Temperature Range	-50°C to +250°C
Cable Length	1m, 2m, 5m
Cable Type	7/0.2mm Teflon® insulated 2 core twisted leads
Termination Type	2-wire tails
Accessories	1/8" BSPP brass male compression fitting with locknut included