



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

IEC Screw In Thermocouple Temperature Probes - Type K - 1/4" BSPP - 100mm Length x 3mm Diameter



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Screw in type K temperature sensors with a measuring range of -75 to +250°C and 1 metre PFA flat pair lead wire.

For direct mounting into vessels, pipes, tanks, or ducting for liquid or gas temperature measurement.

Typical Applications

Industrial Process Control

Pipes or Ducts

Vessels & Tanks

Specifications

Specifications

Product Code	XE-7333-001
General Description	Screw in Type K thermocouple
Thermocouple Type	K IEC
Cable Length	1M
Cable Type	PFA
Termination Type	Stripped Wires
Sheath Length (mm)	100mm
Sheath Diameter (mm)	3mm
Process Connection options	1/4"BSPP
Max. Temperature	250°C
Min. Temperature	-75°C
Standards Met	IEC
Number of Wires	2