



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

IEC Fabricated Tube Thermocouples with Crimp Seal and Lead Wire - Type T- 6mm Sheath Diameter - 100mm Sheath Length



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Typical Applications

Industrial Process Control

Pipes or Ducts

Ovens and heating

General purpose

Specifications**Specifications**

| | |
|-----------------------------|---|
| Product Code | XE-7933-001 |
| General Description | Fabricated IEC Thermocouple with Rigid Stainless Steel Sheath with 1 metre PFA lead |
| Thermocouple Type | T |
| Cable Length | 1M |
| Cable Type | PFA |
| Termination Type | Stripped Wires |
| Sheath Length (mm) | 100mm |
| Sheath Diameter (mm) | 6mm |
| Max. Temperature | 250°C |
| Min. Temperature | -75°C |
| Tolerance | Class 1 |
| Standards Met | IEC |
| Number of Wires | 2 |
| Series | Fabricated Tube Thermocouple |