



LABFACILITY
TEMPERATURE & PROCESS TECHNOLOGY

IEC Fabricated Tube Thermocouples with Crimp Seal and Lead Wire - Type T- 4.5mm Sheath Diameter - 150mm Sheath Length



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Typical Applications

Industrial Process Control

Pipes or Ducts

Ovens and heating

General purpose

Specifications**Specifications**

Product Code	XE-7930-001
General Description	Fabricated IEC Thermocouple with Rigid Stainless Steel Sheath with 1 metre PFA lead
Thermocouple Type	T
Cable Length	1M
Cable Type	PFA
Termination Type	Stripped Wires
Sheath Length (mm)	150mm
Sheath Diameter (mm)	4.5mm
Max. Temperature	250°C
Min. Temperature	-75°C
Tolerance	Class 1
Standards Met	IEC
Number of Wires	2
Series	Fabricated Tube Thermocouple