



## Surface Mount RTD Sensor with Miniature 3-PIN plug - 1M Length



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

A bolt-on RTD Temperature Sensor with an aluminium housing for fast thermal response and tab with 3mm mounting hole for easy mounting, complete with a PFA lead and tails.

The sensor can be cemented, glued, screwed, or bolted on most flat surfaces.

It has a Pt100 Class A sensor as standard with optional Pt1000 sensor and lead lengths to suit, as well as non-standard connectors options, such as Lemo plug to suit your application.

### Specifications

## Specifications

<b>Product Code</b>	XE-7703-001
<b>General Description</b>	RTD-830 style - Stripped 75mm tails as standard
<b>Termination Type</b>	Miniature 3-PIN Plug
<b>Max. Temperature</b>	230°C
<b>Min. Temperature</b>	-50°C
<b>Tolerance</b>	Class A
<b>Length</b>	1M Lead
<b>Number of Wires</b>	4 wire or 3 wire with connector
<b>Series</b>	Bolt on RTD