




LABFACILITY
TEMPERATURE & PROCESS TECHNOLOGY

Surface Mount RTD Sensor without plug - 1M Length



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

A bolt-on RTD Temperature Sensor with an aluminium housing for fast thermal response and tab with 3mm mounting hole for easy mounting, complete with a PFA lead and tails.

The sensor can be cemented, glued, screwed, or bolted on most flat surfaces.

It has a Pt100 Class A sensor as standard with optional Pt1000 sensor and lead lengths to suit, as well as non-standard connectors options, such as Lemo plug to suit your application.

Specifications

Specifications

Product Code	XE-7700-001
General Description	RTD-830 style - Stripped 75mm tails as standard
Termination Type	Bare Tails
Max. Temperature	230°C
Min. Temperature	-50°C
Tolerance	Class A
Length	1M Lead
Number of Wires	4 wire
Series	Bolt on RTD