


## Surface Mount Type K Thermocouple Sensor



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Surface Mount Type K Thermocouple.**

A bolt-on type thermocouple with an aluminium housing for fast thermal response and tab with 3mm mounting hole for easy mounting. It is a Class 1 accuracy type K thermocouple as standard complete with a PFA lead with tails or mini plug. The sensor can be cemented, glued, screwed, or bolted on most flat surfaces.

**Specifications****Specifications**

<b>Sensor Type</b>	Thermocouple Type K
<b>Fitting</b>	3mm Mounting Hole
<b>Cable Type</b>	PFA
<b>Termination Type</b>	Bare Tails or Miniature Plug Options
<b>Max. Temperature</b>	230°C
<b>Min. Temperature</b>	-50°C
<b>Tolerance</b>	Class 1
<b>Length</b>	1M, 2M or 5M Lead
<b>Important Information</b>	Optional Type J or T version available on request