

## Pipe Clamp Surface Pt100 RTD Jubilee Style with PFA Lead - 13-20mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

A surface temperature sensor with a measuring range of -50 to +250°C with a Pt100 class A 4 wire resistance thermometer detector (RTD).

Available in four sizes with adjustable jubilee clip to fit pipe diameters of 13-20mm, 20-32mm, 30-45mm and 40-55mm.

Designed to be easily installed for robust clamp contact with pipes or round surfaces with various diameters.

With a 2 metres PFA insulated flexible lead wire and stripped tails.

Pt100 RTD class A 4 wire resistance thermometer as per IEC 751

Sensor temperature rating -50 to +250°C

Adjustable to suit pipe diameters 13-20mm, 22-30mm, 30-45mm and 40-55mm

### **Specifications**

## **Specifications**

<b>Product Code</b>	XE-7790-001
<b>Sensor Type</b>	Pt100
<b>Cable Length</b>	2M
<b>Cable Type</b>	PFA
<b>Cable Termination</b>	Stripped wires
<b>Max. Temperature</b>	250°C
<b>Min. Temperature</b>	-50°C
<b>Tolerance</b>	Class A
<b>Number of Wires</b>	4 wire